

# Cleaning Your Nebuliser

Regular care of the nebuliser system is very important. To maintain the proper operation of your nebuliser, it is recommended that you clean and store it according to the manufacturer's instructions. It is important to note that different nebuliser systems may require different cleaning procedures.

**Clean your nebuliser and its pieces after each use.**



1 Wash the container and mask or mouthpiece with dish soap and hot water.

When disinfecting (every other day), boil the nebuliser parts for approximately 10 minutes except for PVC parts (some masks and tubes). For PVC parts, alternatively, use commercially available disinfectant.



2 Shake off the excess water and let the parts air dry.



3 You can also attach the nebuliser pieces to the machine and turn it on to dry the nebuliser quickly.



4 Make sure all pieces are completely dry before storing them away.

Always clean and care for your nebuliser so that it can take care of you.